



LEMO 10 Pin-M Connector 1B - FGG (frontal view)



- 1--- GND VUELTA
- 2--- FN
- 3--- FL
- 4--- SUM
- 5--- +12V (white)
- 6--- GND
- 7--- -12V
- 8--- V LASER
- 9--- BIAS
- 10--- V TAPPING

LEMO 10 Pin-M Connector 1B - FGG (back view)



- 1--- GND VUELTA
- 2--- FN
- 3--- FL
- 4--- SUM
- 5--- +12V (white)
- 6--- GND
- 7--- -12V
- 8--- V LASER
- 9--- BIAS
- 10--- V TAPPING