



LEMO Connector FGG-1B 308 (8 Pin-P)



Frontal view

- 1-----+Y
- 2-----X
- 3-----+Z (outside external piezo)
- 4-----Z (inside internal piezo)
- 5-----+X (white)
- 6-----GND
- 7-----Y
- 8-----Z (inside external piezo)

LEMO Connector FGG-1B 308 (8 Pin-P)



Back view

- 1-----+Y
- 2-----X
- 3-----+Z (outside external piezo)
- 4-----Z (inside internal piezo)
- 5-----+X (white)
- 6-----GND
- 7-----Y
- 8-----Z (inside external piezo)