


LEMO 10 Pin-P/ Connector FGG-1B Male (frontal view)


AFM / STM



- 1--- GND power
- 2--- FN / I (current)
- 3--- FL / Reserved
- 4--- SUM / Reserved
- 5--- +12V (white)
- 6--- GND
- 7--- -12V
- 8--- V Laser / Reserved
- 9--- BIAS
- 10--- V Dynamic / Reserved

LEMO 10 Pin-P/ Connector FGG-1B Male (back view)

AFM / STM



- 1--- GND power
- 2--- FN / I (current)
- 3--- FL / Reserved
- 4--- SUM / Reserved
- 5--- +12V (white)
- 6--- GND
- 7--- -12V
- 8--- V Laser / Reserved
- 9--- BIAS
- 10--- V Dynamic / Reserved